

DRE Vista Flex LED

Examination Light

Flexible LED examination light offers variety of mounting and illumination options

The DRE Vista Flex is a versatile LED examination light with maximum positioning capabilities, making it ideal for Dermatology, Plastic Surgery, ENT, OB/GYN, Ophthalmology, Outpatient Surgery, Family Practice and more. The dual-segmented gooseneck arm and universal mounting clip allow clinicians to position the light precisely where they need it.

Complete with adjustable spot, user-friendly control panel and an extended length power cable.

Features

- ✓ Bright LED examination light ideal for a wide range of medical specialties.
- ✓ Cost-saving, energy-efficient LED bulb produces up to 50,000 hours of illumination.
- ✓ Available in a variety of mounting options including a mobile stand, wall mount or table mount.
- ✓ Universal clip allows clinicians to easily attach the light to any table or procedure chair.
- ✓ Flexible gooseneck light arm features two segments for maximum positioning capabilities.
- ✓ Light produces 70,000 Lux at an eight inch working distance.
- ✓ Bright white color temperature accurately illuminates tissue.
- ✓ Adjustable spot ranges between 0.4 to 8.6" (depending on distance from workspace).
- ✓ Control panel features intuitive controls and two handles for easily maneuvering the lighting unit.
- ✓ Streamlined design for easy cleaning and sterilization.
- ✓ Power cable measures nearly 10 feet.



SPECIFICATIONS



LED: 1 White 3W LED Bulb

Lifetime: 50,000 hours

Color Temperature:
5,300K

Adjustable Spot Diameter
(working distance 7.87"):
0.39 - 1.7"

Illuminance
(working distance 7.87"):
70,000 lux



Gooseneck Length: 39.37"

Stand Height: 27.56"

Base Diameter: 19.69"

Package Measurement:
33.86 x 24.02 x 6.3"



Gross Weight: 6kg (13.22lb)

Net Weight: 3.5kg (7.72lb)



Voltage:
100-240V AC ~ 50/60Hz

Power: 10 VA

Power Cable:
#18 AWG, 118.11"

Fuse:
T2AI/250V

1393-01-2016-12