

DRE Meridian

Power Exam Table

Designed for patient comfort and convenience

The DRE Meridian Exam Table is an economical option for facilities seeking a simply designed, convenient exam table. The table comes with a pull-out foot rest and extendable stirrups, making it ideal for a variety of office-based medical exams. With electronic height adjustment, adjustable foot stirrups, and convenient paper roll holder, the Meridian Exam Table can be easily prepared between patients in any busy practice.

Features

- ✓ Ergonomically designed seat section and back section are designed for patient comfort.
- ✓ Seat section features pelvic tilt for enhanced positioning capabilities.
- ✓ Convenient and flexible operation with quiet motor.
- ✓ Easy-to-clean, smooth surface.
- ✓ Foot stirrups angle can be positioned for each patient: 15°, 30°, 45°
- ✓ Extendable, padded foot rest.
- ✓ Use convenient handle to manually adjust back and seat rest.
- ✓ Electronic height adjustment via hand control.
- ✓ Stainless steel basin slides out from the base of the table.
- ✓ Convenient paper roll holder located under the back section.
- ✓ Height adjustment range of 19.5 - 39 inches.



1562-2017-09-12

SPECIFICATIONS



Maximum load:
380 lb



Dimensions:
55 in ± 0.8 in x 28 in ± 0.8 in

Height:
19.5 to 39 in

Slide length of foot rest:
11 in ± 0.8 in

Slide length of pull-out basin:
17 in ± 0.8 in

Table height adjustment:
19 in ± 0.4 in



Back section angle:
0° ~ 70° ± 5°

Seat section angle:
0° ~ 7° ± 1°

Foot rest angle:
45° ± 5°