



MacIntosh Blades

English Profile


Laryngoscopes

With Fiber Optic
Illumination



MACINTOSH - ENGLISH (PENLON™ TYPE)

The Original English profile with a soft, gentle curvature and a flange starting at the lock of blade and continuing to the tip as it narrows. Hand crafted in soft, matte finish surgical stainless steel to virtually eliminate reflection and glare.



All above blades are made from 303/304 surgical stainless steel and are provided with 4mm, replaceable fiber optic bundles. GreenLine® blades and handles are fully compliant with ISO standard 7376-3 and interchange perfectly with companies conforming to the "Green System". (i.e. Anesthesia Associates, AMS, Heine, Propper, Rusch and Welch Allyn)

CAT#	SIZE	DESCRIPTION	A	B
5-5232-01	1	Infant	92mm	22mm
5-5232-02	2	Child	100mm	21mm
5-5232-03	3	Medium Adult	130mm	22mm
5-5243-03	3	Left Hand Adult	130mm	22mm
5-5232-04	4	Large Adult	155mm	25mm
5-5232-05	5	Extra Lg. Adult	175mm	22mm